
Combi

Combi Multifunctional Compact Thermometer User Manual

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Thank you for purchasing the Combi Multifunctional Compact Thermometer. Please read this manual carefully before use and use accordingly. Please keep this manual for future reference. If you transfer this product to others, please also attach this manual to the product.

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Index

Before use.....	1
To use it safely.....	1
Names of parts and sales package.....	2
Operation.....	2
Method of battery replacement.....	3
Troubleshooting.....	4
Maintenance and storage.....	4
Materials and specification.....	4

Before use

- This product is able to measure the forehead temperature in 1 second.
- This product can measure surface temperature. It helps monitoring if the object temperature is suitable for the baby, for example milk or bath water.
- The protection cap is used to keep the thermometer clean.
- Please do not use for other than the above purpose.
- Pull out the insulated strip from the battery compartment before use.

Age reference: From newborn

To use it safely

• When using

CAUTION

- Do not use this product other than for the purpose mentioned.
- Never leave child unattended while using this product.
- Do not place the product near fire.
- Do not drop or bang the product as it may cause damage.
- Check for damage before every use. Stop using this product immediately, if any crack or damage is found.
- Keep thermometer out of child's reach.
- Do not immerse the thermometer into any liquid, otherwise it may be damaged.
- Do not try to remodel, disassemble, or repair the product.
- Stop using the thermometer if it produces a strange odor, makes a strange sound or smoke appears.
- Do not place the thermometer close to any object that emits strong static electricity or electromagnetic wave such as television, microwave oven, mobile phone etc. Otherwise, product may be damaged.
- Keep battery out of child's reach. If child swallows it, contact doctor immediately.
- Battery must be replaced by adult.
- Misuse of battery may cause it to leak, overheat or explode. People may get injured, product may be broken and an accident may occur. Therefore, please note the following cautions.
- Read battery warning and caution, and use accordingly.
 - Remove battery when not in use for extended periods.
 - Wash the clothes immediately if liquid from battery leaks.
 - In the event battery liquid comes into contact with eyes or skin, rinse the affected area thoroughly with water, and see doctor immediately.
- The product is intended for household use only.
- Do not measure body temperature within 30 minutes after meals, exercise or bath.
- When the environmental conditions change, for example from indoor to outdoor, please stay in the new environment for about 30 minutes before measuring the new temperature.


• Storage

CAUTION

- Always put a protection cap on the probe for storage.
- If the cap gets wet, please dry it before putting on the probe. Otherwise, the thermometer may become damaged.
- Do not leave the product in direct sunlight, in the car, in a dusty area or near water. Otherwise, it may become damaged.

Names of parts and sales package

• Sales packaging:

- Combi Non-contact Forehead thermometer
- Protection Cap
- User manual

Operation


Use of the thermometer

1 Forehead temperature

The thermometer has been designed for practical use. It is not meant to replace a visit to the doctor. Please also remember to compare the measured result to your regular body temperature. Please consult with doctor if you have health concerns.

Temperature measuring operation:

Press the  button to power on the device. Forehead mode is the default mode. You can see the  icon on the screen. In this mode, you can press the "Scan" button and hear one beep sound to get the body temperature by forehead measurement.

1.1 Power on:

Press the  button to switch the device on. The device performs a brief self test: all elements of the display are shown on the screen.

1.2 Measuring body temperature on the forehead:

Hold the thermometer 1 to 3 cm away from the central forehead and press the "Scan" button. While scanning, you will hear a beep sound, which indicates you the newest measurement is taking place. The time consuming for measurement might be 1 second.

※ If body temperature is over 37.5°C, fever alarm will be heard with a long beep sound following with 3 short beeps.

1.3 For a second reading, wait until the icon stops flashing on the display with 2 beep sounds. Then, it is ready for the next measurement.

CAUTION

- Forehead temperature is displayed in oral mode. This mode converts the forehead temperature to display its "oral-equivalent" value.
- Before the measurement, please stay in a stable environment for 5 minutes and avoid exercise, or bath for 30mins.
- Remember to keep the forehead area clean and away from sweat, cosmetics and do not scan on the scar position while taking temperature.

WARNING

- Keep the protection cap away from children while in use.

2. Surface temperature

The surface mode shows the actual and unadjusted surface temperature which is different from the body temperature. It can help you monitor if the object temperature is suitable for baby, for example baby's milk or bath water.

Temperature measuring operation:

After power on, press and hold the  button, and simultaneously press the "SCAN" button one time for "Infrared thermometer" mode to see  icon on your LCD display. In this mode, you can get the target surface temperature.

※If you want to change it back to Forehead Mode, repeat the above steps.


2.1 Measuring surface temperature:

- When you press the "Scan" button, you will get the real time temperature immediately. If you press and hold the "Scan" button, the reading of measurement will be continuously updated.
- Applications include temperature measurements for water, milk, cloth, skin or other objects

Note: This mode shows the actual and unadjusted surface temperature which is different from the body temperature.

Other functions:


Room temperature

Suitable ambient temperature is important for the baby. The thermometer always helps you recognize the room temperature. Temperature measuring operation:

The room temperature will be automatically shown on the screen with  icon, after left it idle for about 1 minute from Forehead mode or Surface mode.

Memory locations

There are total 25 sets of measurement condition for forehead measurements.

When power on, Press the  button to see the temperature records with  icon.

Last reading

After getting a new temperature reading, the last reading will be shown on the screen (at top right corner) with the  icon.


Power Off

The temperature automatically turns to room temperature if left idle for more than 1 minute.

At temperature measuring mode, press and hold  button for 4 seconds, the device would power off.

Method of battery replacement

1. Insert a pointed object, such as a toothpick, into battery cover pin hole. At the same time, use thumb to push battery cover out.
2. Hold the device and flip the battery out with a small screwdriver.
3. Insert the new battery (CR2032) under the metal hook on the side, and press the battery down until you hear a click.
4. Slide the battery cover back until it is locked.

WARNING

- Keep battery out of child's reach. If child swallows it, contact doctor immediately.

CAUTION

- Inside battery with the positive (+) side is up facing while the negative (-) is down facing.
- Dispose the battery at appropriate collection sites according to national or local regulations.

Troubleshooting

Error Message	Problem	Solution
Er	Error 5~9, the system is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
Er1	Measure before device stabilization.	Wait until all the icons stop flashing.
Er2	The device shows a rapid ambient temperature change.	Allow the thermometer to rest in a room for at least 30 minutes at room temperature: 10°C and 40°C.
Er3	The ambient temperature is not within the range between 10°C and 40°C.	
Hi	(1) In Forehead mode: Temperature taken is higher than +42.2°C (2) In Surface mode: Temperature taken is higher than +80°C	Please select the target within specifications. If a malfunction still exists, please contact the nearest retailer.
Lo	(1) In Forehead mode: Temperature taken is lower than +34°C (2) In Surface mode: Temperature taken is lower than -22°C	
	Device cannot be powered on to the ready stage.	Change with a new battery.

Maintenance and storage

- The probe lens is the most delicate part of the thermometer. Use with care when cleaning the probe lens to avoid damage.
- Use the 70% alcohol swab or the soft cotton moistened with the 70% alcohol to clean the probe lens.
- Allow the probe to fully dry for at least 1 minute.
- Keep the unit dry and away from any liquids and direct sunlight.
- Storage temperature range: It should be stored at room temperature between -20~+50°C, RH 85%
- The probe should not be submerged into liquids.
- Holding the thermometer too long may cause a higher ambient temperature reading of the probe. This could make the body temperature measurement lower than usual.

Materials and Specifications

Materials

Body of thermometer : ABS
Protection cap : ABS

Battery

1pc lithium cell battery (CR2032)

Specifications

Temperature measurement range:

- Forehead mode: 34~42.2°C
- Surface mode: -22~80°C

Operating temperature range: 10~40°C

Storage temperature range: It should be stored at room temperature between -20~+50°C, RH ≤ 85%

Transportation temperature shall be less than 70°C, RH ≤ 95%, Comply with ASTM E1965-98, IEC/EN60601-1-2(EMC), IEC/EN60601-1(Safety) standards.

Accuracy: • Forehead mode: ±0.2°C (in measured range at 36.0~39.0°C), ±0.3°C (in measured range at 34.0~35.9°C and 39.1~42.2°C).

- Surface mode: ±4% of reading or ±2°C whichever is greater.

Battery life is estimated at 3,000 readings. Enable normal use: 1 year.

This thermometer converts the forehead temperature to display its "oral equivalent."

CAUTION

- There is no gender and age limitation for using the infrared thermometer.
- This is not an AP or APG product.

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